

F

B'SYS GmbH

# CHO-Tet hERG Cell Line

Specification Sheet

© B'SYS GmbH

---

## TABLE OF CONTENTS

<b>1</b>	<b>BACKGROUND</b>	<b>3</b>
1.1	Drug-induced QT Prolongation	3
1.2	Regulatory Issues	3
1.3	B'SYS' CHO-Tet hERG Cells	3
<b>2</b>	<b>PRODUCT SHIPMENT</b>	<b>4</b>
2.1	Product Format	4
2.2	Mycoplasma Certificate	4
<b>3</b>	<b>VALIDATION OF CHO-TET HERG CELLS</b>	<b>5</b>
3.1	Electrophysiology	5
3.2	Tail Current Activation	5
3.3	Inward Rectifying	6
3.4	Positive Control E-4031	6
3.5	Patch-clamp Success Rates	6
<b>4</b>	<b>CELL CULTURE CONDITIONS</b>	<b>7</b>
4.1	General	7
4.2	Recommended Complete Medium	7
4.3	Antibiotics	7
4.4	Thawing Cells	7
4.5	Splitting Cells	7
4.6	Freezing Medium	8
4.7	Freezing Cells	8
4.8	Induction of channel expression	8
4.9	Stability of CHO-Tet hERG cells	8
<b>5</b>	<b>HERG SEQUENCE</b>	<b>9</b>
5.1	Accession Number U04270	9
<b>6</b>	<b>CONTACT INFORMATION</b>	<b>10</b>
6.1	Contact Address for Technical Support & Ordering Information	10

---

## 1 BACKGROUND

### 1.1 Drug-induced QT Prolongation

The QT interval as measured in the electrocardiogram (ECG) is determined by the duration of the cardiac action potential. The rapid delayed rectifier current ( $I_{Kr}$ ) is important for cardiac action potential repolarization. Suppression of  $I_{Kr}$  function by adverse drug effects (e.g. terfenadine, cisapride, astemizole etc.) can induce a prolongation of the QT interval carrying elevated risk of life-threatening cardiac arrhythmias.

### 1.2 Regulatory Issues

$I_{Kr}$  is mediated by the potassium channel hERG (human ether-à-go-go related gene). This ion channel is involved in both, congenital and acquired long QT syndrome. According to the ICH S7A and S7B guidelines the analysis of molecules for inhibition of hERG currents allows the assessment of the potential for QT prolongation.

### 1.3 B'SYS' CHO-Tet hERG Cells

As a new test system B'SYS has designed a CHO-Tet hERG cell line with inducible expression of hERG channels for extra long culturing without loss of channel expression. As special feature channel expression can be easily fine-tuned for single-channel recordings. The hERG cDNA was cloned and transfected into CHO cells and then the functional properties of the hERG channels validated by means of the patch-clamp technique.

## 2 PRODUCT SHIPMENT

### 2.1 Product Format

CHO-Tet cells stably transfected with recombinant hERG potassium channel:

- 1 x 0.5 mL aliquots of frozen cells at 2.7 E+06 cells/mL
- Cells are frozen in complete medium with 10% DMSO
- Cells are frozen at passage number 15

### 2.2 Mycoplasma Certificate

B'SYS periodically tests cells for presence of mycoplasma by means of highly sensitive PCR based assays. All delivered cells are free of mycoplasma.

### 3 VALIDATION OF CHO-TET HERG CELLS

#### 3.1 Electrophysiology

hERG currents were measured by means of the patch-clamp technique in the whole-cell configuration. The bath solution contained (in mM) NaCl 137, KCl 4, CaCl<sub>2</sub> 1.8, MgCl<sub>2</sub> 1, D-glucose 10, HEPES 10, pH (NaOH) 7.40. The pipette solution consisted of (in mM) KCl 130, MgCl<sub>2</sub> 1, MgATP 5, HEPES 10, EGTA 5, pH (KOH) 7.20. After formation of a GΩ seal between the patch electrodes and individual hERG stably transfected CHO-Tet cells, the cell membrane across the pipette tip was ruptured to assure electrical access to the cell interior. All solutions applied to cells were continuously perfused and maintained at room temperature. As soon as a stable seal could be established hERG currents were measured upon depolarization of the cell membrane from a holding potential of -80 mV to +60 mV in 10 mV increments of 2 s duration (Fig. 1). hERG tail currents were elicited upon partial repolarization to -40 mV for 3 s. The voltage pulses were run at intervals of 10 s.




Fig. 1: Representative hERG outward currents recorded upon depolarization of the CHO membrane from -80 mV to +60 mV in 10 mV increments (2 s). Subsequent partial repolarization to -40 mV (3 s) elicited hERG tail currents.

#### 3.2 Tail Current Activation

The hERG channel is known as an inwardly rectifying potassium channel displaying very slow activation/deactivation kinetics, but is inactivated very rapidly. To verify the functional properties of B'SYS's cloned CHO-Tet hERG cells, both the activation plot and IV curve were constructed. The activation plot of hERG tail currents was best fitted with a Boltzmann function with half-maximal activation at **-9.2 mV** and a slope of **8.6** (Fig. 2, n=9 cells).




Fig. 2: Activation plot of hERG tail currents fitted to Boltzmann function with V<sub>50</sub> = -9.2 mV, k = 8.6 and max = 1623 pA (mean of 9 cells).




Fig. 3: IV curve of hERG activation currents. Threshold for activation was between -40 mV and -30 mV and maximal current amplitude was measured at 0 mV (mean of 9 cells).

### 3.3 Inward Rectifying

Smallest hERG activation currents could be elicited upon depolarization of the cell membrane to around -30 mV (Fig. 3). Further depolarization gradually increased the outward currents until a maximum was reached at 0 mV. Depolarization of the membrane to positive potentials resulted in a decrease of current amplitudes due to the inward rectifying properties of the hERG channel.

### 3.4 Positive Control E-4031

E-4031, which is known as a selective  $I_{Kr}$  blocker, was used as positive control to validate the pharmacological profile of the CHO-Tet hERG cells. The effect of E-4031 was first measured at 30 nM, which blocked hERG tail currents almost completely ( $8.9 \pm 3.6$  % relative tail current amplitude, mean  $\pm$  SEM of 3 cells). When the cells were exposed to 10, 3 and 1 nM E-4031, the steady-state level relative tail current amplitudes revealed values of  $27.9 \pm 2.8\%$ ,  $77.8 \pm 3.1\%$  and  $99.4 \pm 2.2\%$ , respectively (mean  $\pm$  SEM of 3 cells). The inhibition curve (Fig. 4) was best fitted with an  $IC_{50}$  value of **6.01 nM**.




Fig. 4: hERG tail current inhibition curve of E-4031.

### 3.5 Patch-clamp Success Rates

The patch-clamp properties of the CHO-Tet hERG cell line were elucidated at typical working passage numbers (passage 12-44). A total of 14 cells were analyzed. Success for establishment of on-cell configuration was defined as follows:  $> 1$  G $\Omega$ . The whole-cell configuration was not accepted if the membrane resistance was below 500 M $\Omega$ . A successful recording had to be free of rundown effects and variations in series resistance.

- On-cell successful: **91%** (n=14)
- Whole-cell successful: **86%** (n=14)
- Recording (20-30 min) successful: **76%** (n=14)

## 4 CELL CULTURE CONDITIONS

### 4.1 General

CHO-Tet hERG cells are incubated at 37°C in a humidified atmosphere with 5% CO<sub>2</sub> (rel. humidity > 95%). The cells are continuously maintained and passaged in sterile culture flasks containing F12 (HAM) medium supplemented with 9% fetal bovine serum, 0.9% Penicillin/Streptomycin solution, 100 µg/mL Hygromycin and 15 µg/mL Blastidicin. The CHO-Tet hERG cells are passaged at a confluence of about 80%. For electrophysiological measurements the cells are seeded onto e.g. 35 mm sterile culture dishes containing complete medium.

- All solutions and equipment coming in contact with the cells must be sterile.
- Use proper sterile technique and work in a laminar flow hood.
- Be sure to have frozen cell stocks at hand before starting experiments.
- Cells should be split every 2-3 days at 70% - 80% confluency at 1:3 to 1:5 ratio.

### 4.2 Recommended Complete Medium

- F12 (HAM) with GlutaMAX I
- 9% FBS
- 0.9% Penicillin/Streptomycin

### 4.3 Antibiotics

- CHO-Tet hERG clones were selected under 500 µg/mL Hygromycin and 100 µg/mL Blastidicin pressure.
- To cultivate CHO-Tet hERG cells, a reduced antibiotic pressure (see 4.1.) must be used.
- To separate CHO-Tet hERG cells from untransfected cells, use the selection pressure.

Remark: The permanent application of antibiotics has no effect on current density.

### 4.4 Thawing Cells

- Remove vial of cells from liquid nitrogen and thaw quickly at 37°C.
- Decontaminate outside of vial with 70% ethanol.
- Transfer cells to a T-25 culture flask containing 5 mL complete medium.
- Incubate cells at 37°C for 4-6 h to allow the cells to attach to the bottom of flask.
- Aspirate off the medium and replace with 5 mL complete medium & antibiotics.
- Antibiotics: 100 µg/mL Hygromycin and 15 µg/mL Blastidicin.
- Incubate cells and check them daily until 70% - 80% confluency is reached.

### 4.5 Splitting Cells

- When cells are 70% - 80% confluent remove medium.
- Wash cells once with 1x PBS to remove excess medium.
- Add 1x Trypsin/EDTA and incubate 30 s at room temperature.
- Remove Trypsin/EDTA quickly and incubate cells for 2 min at 37°C.
- Detach cells, add complete medium and pipet up and down to break clumps of cells.
- Passage cells into new flask with complete medium and antibiotics at 1:3 to 1:5 ratio.
- Use remaining suspension for counting the cells.

#### 4.6 Freezing Medium

- Mix 0.9 mL fresh complete medium and 0.1 mL DMSO for every 1 mL freezing medium.
- Sterilize freezing medium by means of appropriate micro filter (0.1 µm – 0.2 µm).

#### 4.7 Freezing Cells

- Prepare fresh freezing medium and keep it on ice.
- Cells should have 80% - 90% confluency prior to freezing.
- Remove the complete medium.
- Wash cells once with 1x PBS to remove excess medium.
- Add 1x Trypsin/EDTA and incubate 30 s at room temperature.
- Remove Trypsin/EDTA quickly and incubate cells for 2 min at 37°C.
- Detach cells, add complete medium and pipet up and down to break clumps of cells.
- Pellet cells with centrifuge and carefully aspirate off medium.
- Resuspend cells at a density of approximately 2.0 E+06 cells per mL with fresh freezing medium.
- Aliquot 0.5 mL of cell suspension into each cryovial.
- Overnight incubate cells in a styropor box at –80°C.
- The next morning transfer cryovial in liquid nitrogen tank for long-term storage.

#### 4.8 Induction of channel expression

CHO-Tet hERG cells are activated by adding 1 µg/mL Tetracycline (final concentration) to the culture medium approximately 24 h prior to hERG experiments.

#### 4.9 Stability of CHO-Tet hERG cells

CHO-Tet hERG cells stably express functionally active hERG potassium channels over 100 passages. Under recommended cell culture conditions no variation in current density was observed over 100 cell splitting cycles.



## 5 HERG SEQUENCE

### 5.1 Accession Number U04270

Cloned cDNA sequence of hERG channel was error-free and identical with U04270 sequence:

```
ATGCCGGTGCAGGGGGCCACGTCGCGCCGAGAACACCTTCCTGGACACCATCATCCGCAAGTTTGAGGGCCAG
AGCCGTAAGTTCATCATCGCCAACGCTCGGGTGGAGAACTGCGCCGTATCTACTGCAACGACGGCTTCTGCGAG
CTGTGCGGCTACTCGCGGGCCGAGGTGATGCAGCGACCCTGCACCTGCGACTTCTGACGCGGGCCGCGCACGCAG
CGCCGCGCTGCCGCGCAGATCGCGCAGGCACTGCTGGGCGCCGAGGAGCGCAAAGTGGAAATCGCCTTCTACCGG
AAAGATGGGAGCTGCTTCCATGCTGCTGGTGGATGTGGTGCCTGTAAGAACGAGGATGGGGCTGTATCATGTTT
ATCCTCAATTTGAGGTGGTGTGGAGAAGGACATGGTGGGGTCCCCGGCTCATGACACCAACCACCGGGGCCCC
CCCACAGCTGGCTGGCCCCAGGCCGCGCCAAGACCTTCCGCCCTGAAGCTGCCCGCGCTGCTGGCGCTGACGGCC
CGGGAGTCGTGCTGCGTGGGGCGCGGGCGGGCCGCGGGGGCCGCTGGTGGTGGAGCTGAGGACCTG
ACGCCCGCGGCAACCAGCAGCAGTCTGCTGGCCCTGGACGAAGTGGACAGCCATGGACAACCACGTCAGGGCTC
GGGCGCGGAGGAGCGGCGTGCCTGGTGGTCCCAGGCTTCCGCCCGCAGCGCGCCCGGCCAGCTCCCATCG
CCCCGGGCGCACAGCCTCAACCCGACGCTCGGGCTCCAGCTGCAGCCTGGCCCGGACGCGCTCCCGAGAAAGC
TGCGCCAGCGTGCGCCGCGCCTGTCGCGCCGACGACATCGAGGCCATGCGCGCCGGGGTGTGCCCCCGCCACC
CGCCACGCCAGCACCGGGCCATGCACCCACTGCGCAGCGGCTTGTCAACTCCACCTCGGACTCCGACCTCGTG
CGTACCGCACCATTAGCAAGATTTCCCAAATCACCTCAACTTTGTGGACCTCAAGGGCGACCCCTTCTTGGCT
TCGCCACCAGTGACCGTGAGATCATAGCACCTAAGATAAAGGAGCGAACCACAATGTCACTGAGAAGGTCAAC
CAGGTCTGTCCCTGGGCGCCGACGTGCTGCCTGAGTACAAGCTGCAGGCACCGCGCATCCACCGCTGGACCATC
CTGATTACAGCCCTTCAAGCCGTGTGGGACTGGCTCATCCTGCTGCTGGTCACTACACGGCTGTCTTACA
CCCTACTCGGCTGCCTTCCCTGCTGAAGGAGACGGAAGAAGGCCCGCCTGCTACCGAGTGTGGCTACGCTGCCAG
CCGCTGGCTGTGGTGGACCTCATCGTGGACATCATGTTTATTGTGGACATCCTCATCAACTTCCGCACCACCTAC
GTCAATGCCAACGAGGAGGTGGTCAGCCACCCCGCCGATCGCCGTCCACTACTTCAAGGGCTGGTTCCCTCATC
GACATGGTGGCCGCCATCCCTTCGACCTGCTCATCTTCGGCTCTGGCTCTGAGGAGCTGATCGGGCTGCTGAAG
ACTGCGCGGCTGCTGCGGCTGGTGCCTGGCGCGGAAGCTGGATCGCTACTCAGAGTACGGCGCGGCCGTGCTG
TTCTTGCTCATGTGCACCTTTCGCTCATCGCGCACTGGCTAGCCTGCATCTGGTACGCCATCGGCAACATGGAG
CAGCCACACATGGAATCACGCATCGGCTGGCTGCACAACCTGGGCGACCAGATAGGCAAACCCCTACAACAGCAGC
GGCCTGGGCGGCCCCCTCCATCAAGGACAAGTATGTGACGGCGCTTACTTACCTTCAGCAGCCTCACCAGTGTG
GGCTTCGGCAACGCTCTTCCCAACCAACTCAGAGAAGATCTTCTCCATCTGCGTCATGCTCATTGGCTCCCTC
ATGTATGCTAGCATCTTCGGCAACGTTGTCGGCCATCATCCAGCGGCTGTACTCGGGCACAGCCCGCTACCACACA
CAGATGCTGCGGGTGCGGGAGTTCATCCGCTTCCACCAGATCCCAATCCCTTGCGCCAGCGCCTCGAGGAGTAC
TTCCAGCAGCCTGGTCTTACACCAACGGCATCGACATGAACGCGGTGCTGAAGGGCTTCCCTGAGTGCCTGCAG
GCTGACATCTGCCTGCACCTGAACCGCTCACTGCTGCAGCACTGCAAACCTTCCGAGGGGCCACCAAGGGCTGC
CTTCGGGCCCCGTCATGAAGTTCAAGACCACACATGCACCGCCAGGGGACACACTGGTGCATGCTGGGACCTG
CTCACCGCCCTGACTTTCATCTCCCGGGGCTCCATCGAGATCCTGCGGGGCGACGTCGTGCTGGCCATCCTGGGG
AAGAATGACATCTTTGGGGAGCCTTGAACCTGTATGCAAGGCCCTGGCAAGTGAACGGGGATGTGCGGGCCCTC
ACCTACTGTGACCTACACAAGATCCATCGGGACGACCTGCTGGAGGTGCTGGACATGTACCTGAGTCTTCCGAC
CACTTCTGGTCCAGCCTGGAGATCACCTTCAACCTGCGAGATAACCAACATGATCCCGGGCTCCCCGGCAGTACG
GAGTTAGAGGGTGGCTTCACTCGGCAACGCAAGCGAAGTTGCTCTTCCGAGGCGCACGGACAAGGACACGGAG
CAGCCAGGGGAGGTGTGCGCTTGGGGCCGGGCGGGGCGGGGCGAGGCGGAGTAGCCGGGGCCGGCCGGGGGG
CCGTGGGGGGAGAGCCCGTCCAGTGGCCCCCTCAGCCCTGAGAGCAGTGAAGGATGAGGGCCAGGCCGAGCTCC
AGCCCCCTCCGCTGGTGCCTTCTCCAGCCCCAGGCCCCCCGGAGAGCCGCGGGTGGGGAGCCCTGATGGAG
GACTGCGAGAAGAGCAGCGCACTTGAACCCCTGTGAGGCGCTTCTCAGGAGTGTCCAACATTTTTCAGCTTC
TGGGGGGACAGTCCGGGGCCGCACTACAGGAGTCCCTCGATGCCCGCCCCACCCCGAGCTTCAACATC
CCCCCTTCCAGCCCGGGTCCGGGCCCCGGGGCGACGTGGAGAGCAGGCTGGATGCCCTCCAGCGCCAGCTCAAC
AGGCTGGAGACCCGGCTGAGTGCAGACATGGCCACTGCTCCTGCAGCTGCTACAGAGGCAGATGACGCTGGTCCC
CCGCTTACAGTGTGTGACCACCCGGGGCTGGCCCCACTTCCACATCCCCGCTGTTGCCCGTACGCCCCCTC
CCCACCTCACCTTGGACTCGCTTCTCAGGTTTCCAGTTTATGGCGTGTGAGGAGCTGCCCCCGGGGCCCCA
GAGCTTCCCCAAGAAGGCCCCACACGACGCTTCCCTACCGGGCCAGCTGGGGGCCCTCACCTCCAGCCCCCTG
CACAGACACGGCTCGGACCCGGGCACTAG
```

## 6 CONTACT INFORMATION

### 6.1 Contact Address for Technical Support & Ordering Information

- B'SYS GmbH  
Technology Center Witterswil  
Benkenstrasse 254  
4108 Witterswil  
Switzerland

Tel: +41 61 721 77 44

Fax: +41 61 721 77 41

Email: [info@bsys.ch](mailto:info@bsys.ch)

Web: [www.bsys.ch](http://www.bsys.ch)